

**Bactericidal activity of GAMA Healthcare Ltd.
biocide determined using the European Standard
Test method BS EN 1276:1997 against: *Acinetobacter
baumanii* NCIMB 9216.**

July 2006

Author: P. Humphreys
Authorised by: P. Humphreys

Signature: 
Signature: 

Date: 25/7/06
Date: 25/7/06



University of
HUDDERSFIELD

Tests Carried Out By:

University of Huddersfield School of Applied Sciences
Queensgate
Huddersfield
HD1 3DH

Microbiological Tests

Test Method British/European Standard BS EN 1276:1997.
Dilution-neutralisation

Test Procedures Full details of all the test and control procedures
used are given in the Test Method

Disinfectant GAMA Healthcare Ltd biocide
Batch number: N/A
Date of delivery: June 2006
Storage conditions: 20°C – 25°C
Active substances: not specified
Appearance product dilutions: colourless, clear
product solution.

Interfering Substance (Organic Challenge)

1. Simulated clean conditions:
0.3 g l⁻¹ bovine albumin (final
concentration)
2. Simulated dirty conditions:
3.0 g l⁻¹ bovine albumin (final
concentration)

Temperature Ambient (25°C)

Contact Time Tested 5 (± 10 s) minute.

Test Organisms *Acinetobacter baumannii* NCIMB 9216.

Culture Medium Tryptone Soya Agar, Lab M

Incubation Plates were incubated at 25 °C for 48 h

Diluent MRD, Lab M

Neutraliser Neutraliser, containing 60g/l Tween 80, 60g/l
Saponin, 2g/l L-histidine, 2g/l L-cysteine in
MRD.

General Method

A standard suspension of test organisms containing $1.5 - 5.0 \times 10^8$ cells ml^{-1} of bacteria was prepared. 1 ml of interfering substance was pipetted into a Universal bottle, followed by 1 ml of test organism suspension. The mixture was mixed and left for 2 minutes. After 2 minutes 8 ml of GAMA Healthcare Ltd biocide was added and mixed. After a contact time of 5 minutes, a 1 ml sample of the reaction mixture was pipetted into 9 ml of neutraliser and left for 5 minutes. A 1 ml sample was then pipetted into 2 Petri dishes and mixed with 15 ml of culture medium tempered at 47 °C. After setting, the Petri dishes were incubated at 25°C. Colony forming units were counted after 2 days incubation and the fraction of surviving organisms calculated.

Requirements of this standard

The product, when tested as stipulated under simulated clean conditions (0.3 g l^{-1} bovine albumin) or dirty conditions (3 g l^{-1} bovine albumin) under the required test conditions (25°C, 5 minute contact, for the selected reference strain), shall demonstrate at least a 5 \log_{10} reduction in viable counts.

Results¹

Results from the test are summarised in Tables 1 and 2, a full set of results can be found in Table 3.

Test Conditions	Contact Time (minutes)	Log ₁₀ Reduction Achieved
0.3 g l^{-1} (clean)	5	>5 ¹
3.0 g l^{-1} (dirty)	5	>5 ¹

Table 1. Log₁₀ reductions in *A. baumannii* viable counts following a 5 minute exposure to the test material.

Referenced Organism	Starting concentration CFU ml^{-1}	Final concentration CFU ml^{-1} clean 0.3 g l^{-1} Bovine Albumin	Final concentration CFU ml^{-1} dirty 3.0 g l^{-1} Bovine Albumin
<i>Acinetobacter baumani</i> NCIMB 9216	1.9×10^8 (20 ² , 178 ¹ , 23, 20 ²)	Plate count 0, 0. (Actual 6 \log_{10} reduction)	Plate count 0, 0. (Actual 6 \log_{10} reduction)
CFU = colony forming units ¹ viable count of bacterial colonies, 1 ml sample of 10^{-6} bacterial suspension ² viable count of bacterial colonies, 1 ml sample of 10^{-7} bacterial suspension			

Table 2. Reductions in *A. baumannii* viable counts following a 5 minute exposure to the test material.

¹ See Table of results in Appendix 1.

Interpretation of the Results

When tested against *A. baumannii* NCIMB 9216 with a 5 minute contact time a full strength GAMA Healthcare Ltd biocide met the requirements of the Standard at ambient temperature (25°C) under simulated clean and dirty conditions.

Conclusion

According to EN 1276:1997, the batch provided of GAMA Healthcare Ltd biocide possesses bactericidal activity in 5 minutes at ambient temperature (25°C) under clean conditions (0.3g/l bovine albumin) and dirty conditions (3g/l bovine albumin) for referenced strain *A. baumannii* NCIMB 9216.

Signed:

Dr Paul Humphreys
School of Applied Sciences
University of Huddersfield

Appendix 1

Test Organism	VALIDATIONS										Bacterial Test Suspension		Test Procedure Results							
	Bacterial Suspension	Experimental Conditions Validation		Neutraliser Toxicity Control	Dilution Neutralisation Control		Bacterial Test Suspension		Test Procedure Results											
		Clean	Dirty		Clean	Dirty														
A. <i>baumannii</i>	Nv 1.9E+03	Vc 175	206	193	198	Vc 230	195	Vc 239	222	227	223	10-6	202	178	Vc < 15	15	<	15	15	
		Vc 1.9E+02	A	2.0E+02		B	2.1E+02	C	2.3E+02		2.3E+02		10-7	23	20	Na < 1.5E+02	<	1.5E+02		
	Verification of Methodology Pass N is between 1.5E+8 cfu/ml and 5E+8 cfu/ml, N = 1.9E+08 Yes Nv is between 6E+2 cfu/ml and 3E+3 cfu/ml, Nv = 1.9E+03 Yes CLEAN A ≥ 0.05 x Nv when 0.05 x Nv = 9.6E+01 Yes DIRTY A ≥ 0.05 x Nv when 0.05 x Nv = 9.6E+01 Yes B ≥ 0.05 x Nv when 0.05 x Nv = 9.6E+01 Yes CLEAN C ≥ 0.5 x B when 0.5 x B = 1.1E+02 Yes DIRTY C ≥ 0.5 x B when 0.5 x B = 1.1E+02 Yes										Log10 Reductions (cfu/ml)		Clean > 5		Dirty > 5					

Table 3. Testing of *S. typhimurium* (ATCC 14028) the Gamma Health Care Ltd biocide using BS EN 1276:1997.

